



MEDISMART COMPACT AVSU



Medi Smart Delivering innovative and reliable **Medical Engineering**, tailored for the GCC healthcare sector.




EMERGENCY STANDBY MANIFOLD

The MEDISMAST Systems, Emergency Standby Manifold offers a simple standalone unit or may be used as a standby system to a centralized plant or automatic manifold system to provide a source of compressed air or medical gas supply for distribution into a hospital or laboratory application. The unit is supplied as a complete factory tested package, which can be easily tailored to meet specific customer requirements. The system offers a high quality, reliable product, which can be extended to suit individual applications. Systems available as a standard package, which ensures trouble free installation, and reliable operation. Full warranty and technical backup for one year from commissioning, with extended warranties, and extended technical support packages available on request.

Benefits.

- Manufactured to the slandered requirements
- Available as a wall mounted unit with a free-standing option.
- Includes a monitoring and indicating system.
- Forms part of an integrated system which includes manual and automatic options.
- Integrated cylinder rack.
- Easy clean surfaces.
- Includes pressure relief valve, isolation valve and non-return valve.
- Provided as a complete system including all necessary interconnecting pipework and wiring.
- Cylinder connection options to suit Pin Index, Bul Nose, CGA etc.
- Available to suit all medical gases including Oxygen, Nitrous Oxide, Entonox, Medical Air 4 and Surgical Air 7 bar.
- Simple operation ensures trouble free performance.
- Every manifold is fully tested prior to dispatch.





Medismart – Advancing healthcare with trusted, high-quality **Medical Engineering** across the GCC.



- Phone: +971 42 500 178
- Mobile: +971 55 585 2360
- Email: info@medismartglobal.com
- www.medismartglobal.com
- P.O. Box: 231673, Dubai