



Medi Smart Delivering innovative and reliable **Medical Engineering**, tailored for the GCC healthcare sector.




MANIFOLD SYSTEM

Fully Automatic Manifold.



The Automatic Manifold systems is designed to provide aduty and standby gas supply at a constant pressure viatwo banks of gas cylinders with status monitoring Central control panel with a header assembly for medical gas distribution systems in accordance with current HTM standards HTM 02-01, HTM 2022, IS / ISO 7396 - 1 : 2016 Medismart design a unit which requires a minimum of maintenance to ensure the maximum life of the equipment it does require periodic routine maintenance. The Medismart automatic medical gas manifold plant is supplied as a packaged product ready for immediate installation, tested and prepared at the factory for the specific gas indicated and the manifold control panel shall be type tested for electrical safe.

Emergency Manifold.



The Medismart Emergency Standby Manifold system offers a simple standalone unit or may be used as a standby system to a centralized plant or automatic manifold system to provide a source of compressed air or medical gas supply for distribution into a hospital or laboratory application.



Medismart – Advancing healthcare with trusted, high-quality **Medical Engineering** across the GCC.



- Phone: +971 42 500 178
- Mobile: +971 55 585 2360
- Email: info@medismartglobal.com
- www.medismartglobal.com
- P.O. Box: 231673, Dubai